


# SuperEagle®

PUNCTUM PLUG

## The Most Advanced Silicone Punctum Plug Available!



### **Super** patient comfort

- Lowest profile rim
- Soft, low durometer silicone

### **Super** retention

- Patented Tapered Shaft™
- Pivoting wide-flex nose
- Available in three sizes

DESCRIPTION	SIZE	REF
Small	0.4mm to 0.65mm	3130
Medium	0.6mm to 0.85mm	3131
Large	0.8mm and larger	3132

Supplied: 1 ea. sterile silicone plug per single pack, preloaded, on a disposable inserter/dilator.

U.S. Patent Nos. 4,915,684 & 5,334,137 & 5,741,292 & 6,016,806 & 6,027,470 & 6,041,785 & 6,082,362 and other Patents Pending.

CE Marked



8500 Wolf Lake Dr., Suite 110 • P.O. Box 34877 • Memphis, TN 38133 USA  
1-800-222-7584 • (901) 380-7000 • Fax (901) 380-7001  
www.eaglevis.com • info@eaglevis.com